QUICK START GUIDE

GUIDE DE DÉMARRAGE RAPIDE
GUÍA DE INICIO RÁPIDO
GUIDA D'AVVIO RAPIDO
KURZANLEITUNG
BEKNOPTTE HANDLEIDING
SNABBSTARTSGUIDE
PIKAOPAS ALOITUKSEEN
HURTIGSTARTSGUIDE
HURTIGVEILEDNING
GUIA DE INÍCIO RÁPIDO

COREBELT 1.0
COREBELT 3.0
COREBELT 5.0

 SAFETY WARNING

It is strongly recommended to read all the instructions, contra-indications and safety measures carefully before using the device. The complete user manual can be found at www.compex.com

1. DEVICE DESCRIPTION

A On/off button or Pause button during stimulation
B Power and battery LED indicator
C Program button (only available on COREBELT 1.0), to select the program P1 or P2

P1 first LED solid green
P2 first and second LED solid green

D +/- intensity buttons, to increase or decrease the stimulation intensity

E Stimulation LED indicator
F Charger socket, plug in the USB cable to charge the device


A Holster
B Notches to facilitate the positioning of the belt on the back
C Elastic strap to avoid having any loose ends of the belt hanging down
D Metallic pins, electrodes must be positioned over these metal pins
E Self-fastening area to adjust the size of the belt

2. CHARGING THE DEVICE

Make sure the device is charged before using it (green LED).
If the device is not going to be used for over 3 months, ensure that the battery is fully charged.

3. SETTING UP THE BELT

Setup the size of your belt by measuring your waist circumference and seeing which belt number it matches to.

S/M belt	Attachment area	L/XL belt
65 - 67 cm	0	85 - 87 cm
67 - 72 cm	1	87 - 94 cm
72 - 78 cm	2	94 - 100 cm
78 - 84 cm	3	100 - 107 cm
84 - 85 cm	4	107 - 113 cm
-	5	113 - 114 cm

Tips: How to measure the waist circumference?

Stand up straight, place the tape measure at the level of your belly button and wrap it around your waist. Hold the tape measure against your body, not too tight, and take a reading.

Adjust the belt by positioning the internal attachment on belt number found on the previous step.

Then adjust the other side external attachment and bend the belt to get a S-like shape.

Check if the adjustment is okay by placing the center of the back belt on your spinal cord. The notches on the belt facilitate the positioning.

On the front, the plastic belt holster must be centered on your belly button. If not, remove the belt and restart from the beginning to fine tune the belt adjustment.
Once done, do not change these settings unless your body shape changes.

Remove the belt and stick the electrodes grid side down in-line with the outlines. Make sure that the metallic pins are completely covered by the electrodes.

4. DOWNLOADING THE APP

≥ iOS12

≥ Android 8

It is recommended to download and use the app to access all the features of the device.

Pair your device with your smartphone by following the app's guidance. Then select the Objective you want to achieve (COREBELT 3.0 and 5.0) or the Program you want to perform (COREBELT 1.0).

5. STARTING THE SESSION

Before starting a session, remove the plastic cover from the electrodes.
Be sure to use the right electrodes according to your session:
Abs session -> Front electrodes only
Back session -> Back electrodes only
Core session -> Front and Back electrodes

Position the belt around your waist. Firstly, position the back of the belt on your lower back. Then stretch the belt around your waist. The front electrodes should be positioned on each side of your belly button.

Tips: To avoid having any loose ends of the belt hanging down, you can attach the end of the belt to the self-fastening area and slide it under the elastic strap.

Slip the stimulator into the belt holster, make sure it plugs in properly, switch on the device and start your session.
Tips: At the end of the session, you may leave the electrodes on the belt but use the plastic cover to protect them.

N.B. You can also start a session without the use of the app.
For more information please check the complete user manual at www.complex.com.

6. OBJECTIVES/PROGRAMS

	1.0	3.0	5.0
Ab Objectives (only front electrodes used)			
I want to tone my abs		●	●
I want to reduce my waist		●	●
I want to sculpt my abs		●	●
Lower Back Objectives (only back electrodes used)			
I want to reinforce my lower back		●	●
I want to relax my back			●
I want to massage my sore back			●
Core Objectives (front and back electrodes used)			
I want to reinforce my core stability		●	●
I want to work hard to get result 1*		●	●
I want to work hard to get result 2*			●
Programs			
P1 Work my abs (only front electrodes used)	●		
P2 Reinforce my core stability (front and back electrodes used)	●		

* These objectives combine muscle stimulation and active workout at the same time.

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COPYRIGHT © 2022 by DJO, LLC
DJO France S.A.S
DJO France S.A.S · Centre Européen de Frêt, 3 rue de Bethar · 64990 Mouguerre · France
© 2022-05-10 DJO, LLC - 14-0884 Rev D